

OLYMPUS®

Olympus is about life. About photographic innovations that capture precious moments of life. About advanced medical technology that saves lives. About information- and industry-related products that make possible a better living. About adding to the richness and quality of life for everyone. Olympus. Quality products with a **FOCUS ON LIFE**

RESEARCH STEREOMICROSCOPE SYSTEM

SZX

SZX12/SZX9

For Industrial Applications





Today's First Choice: The SZX Research Stereo Microscope System

Setting New Standards In Observation Clarity, System Flexibility And Ease Of Use

Zoom stereo microscopes play an increasingly important and valuable role in the world of modern industry. Wherever they're used, whatever the application, they must offer consistently clear and effective observation of specimens that range widely in both size and material. That's why system flexibility is crucial — along with adequate working distance, wide zoom range, natural image presentation and ergonomic design for easy, comfortable operation.

The SZX Series sets new standards in all these areas, building on the success of the SZH series that preceded it. The key is new technology: optical design advances derived from UIS optics development, and sophisticated system refinements that ensure compatibility with a wide range of image devices. Designed and built with the renowned Olympus commitment to excellence, the SZX Series opens the way to exciting new horizons in zoom stereo microscopy.





CLICK
CLICK STOP

PLAPO 1X PF
OLYMPUS DF PLAPO 1X
JAPAN

SZX12

Zoom ratio: 12.86

Zoom range: 0.7 \times - 9 \times

Max. field diameter (Objective 0.3 \times): 104.8mm


Aperture iris diaphragm: Integrated

SZX-FOF

Stroke: 80mm

Coarse handle stroke per rotation : 36.8mm

Fine handle stroke per rotation : 1.5mm



SZX9

Zoom ratio: 9.05

Zoom range: 0.63 \times - 5.7 \times

Max. field diameter (Objective 0.5 \times): 69.8mm

SZX-FO

Stroke: 80mm

Coarse handle stroke per rotation: 21.2mm



SZX12 And SZX9: World-beating Performances In Image Recording And Operability.



Coarse and fine focusing unit / SZX-FOF

The coaxial coarse/fine focusing knob allows quick, easy focusing without removing hands. The torque of knob rotation can be fine-tuned (a tool is supplied) and can be locked to prevent the image from blurring while observation is in progress. A counter balance is installed in the unit, ensuring that operation remains smooth even when multiple units of photomicrography and video cameras are attached.



Focusing unit / SZX-FO

The torque of knob rotation is adjusted as with the SZX-FOF, and can be locked at a chosen setting.

Highly rigid joint body

All body joints are dovetailed for maximum rigidity. Dual combination light guide (LG-DFI*) or counter balance (SZX-CSP) can be attached.

*merent type of light guide may be offered in each area.

Long 80mm stroke focusing unit

Because the stroke is 80mm long, objectives with different working distances can be used simultaneously

SZX12 zoom body I SZX-ZB12

From 0.7 X - 9X, the SZX12 offers a zoom ratio of 12.86 — currently unmatched anywhere in the world. Equally impressive is the diameter of the observation field (31.43mm) achieved with a 1X objective. The zoom body has an integrated aperture diaphragm which allows increased stopping effects in high magnification observation, helping the SZX12 to deliver outstanding performance in image recording.

SZX9 zoom body / SZX-ZB9

From 0.63 X - 5.7 X, the SZX's zoom ratio of 9.05 puts it at the top of its class and complements the maximum diameter of observation field (34.92mm) achieved with a 1 X objective. This lens performance allows a long working distance with consistently excellent operability. An optional aperture diaphragm unit (SZX-AS) can be attached.



SZX12/SZX9 zoom bodies specifications

	SZX-ZB12	SZX-ZB9
Zoom range	0.7x-9x	0.63X-5.7X
Zoom ratio	12.86	9.05
Numerical Aperture(Ob.1 X)	0.11	0.10
Field number	22	22
Maximum field diameter (Ob.1 X)	31.43mm	34.92mm
Aperture stop unit	incorporated	optional
Click-stop	incorporated	incorporated
Total magnification range	2.1 X-675X	3.15X- 342X

Focusing units specifications

	SZX-FOF	SZX-FO
Coarse movement range per rotation	36.8mm	21.2mm
Fine movement range per rotation	1.5mm	—
Stroke	80mm	80mm
Counter balance	incorporated	optional
Body attachment	dovetail	dovetail
Max. carrying weight	20kg	10kg


The SZX12 — Offering The Exceptional Zoom Ratio Of 1 : 12.86* And Continuous Magnification Change Of 41 X.

The SZX12 is equipped with Olympus' original parfocal objective. When used with the revolving nosepiece SZX-2RE, which allows the attachment of two objectives, the total stepless magnification range of 3.5X - 144X (using 10X eyepiece, DFPLFL0.5XPF, DFPLFL1.6XPF) makes speedy observation simple. Total magnification for the SZX12 is 2.1X - 675X, and for the SZX9 is 3.15X - 342X — both figures being leaders in their respective classes.


* For all zoom stereo microscopes sold worldwide as of August 1987.




Revolving nosepiece



144X



20X




Circuit board


3.5X



800 Lines/mm Resolution: Outstanding Level Of Visibility.



Existing Olympus product



SZX

Pattern on the AF sensor

Combining the DFPLA01×PF for the SZX12 with a SZX-AL20× auxiliary objective sets another new standard, this time in resolution*: 800 lines/mm. Making full use of Olympus expertise in UIS optics technology, the newly designed objective provides high resolution, high contrast and distortion-free images. The lens' apochromatic characteristics eliminate chromatic aberration and ensure excellent color reproduction in all zoom ranges.

* For all zoom stereo microscopes sold worldwide as of August 1997.

Parfocal objective / PF Series

PF series parfocal objectives are used with the SZX12, and are available in the following magnifications: 0.5×, 1.2× and 1.6×. Revolving nosepiece allows easy switching between lenses, and rapid focus by less than a single rotation of the fine focusing knob.

Auxiliary objective / SZX-AL20×

This is a 2.5× lens for attachment at the top of DFPLA01×PF. It enables high magnification observation at a level which conventional zoom stereo microscopes cannot achieve, thus making the SZX ideal for very fine observations, e.g. of semiconductors and magnetic heads.



SZX-AL20×

DFPLA01×PF

Objectives for SZX12



Objective	N.A.	W.D.(mm)	Mag. range
DFPLFL0.3×	0.033	130	0.21 - 2.7
DFPLFL0.45×	0.050	198	0.32 - 4
DFPLFL0.5×PF	0.055	70	0.35 - 4.5
DFPLA01×PF	0.110	74	0.7 - 9
DFPLA01.2×PF	0.132	60	0.84 - 10.8
DFPLFL1.6×PF	0.176	34	1.1 - 14.4
SZX-AL20×	0.275	10/7	10 - 22.5

Objectives for SZX9



Objective	N.A.	W.D.(mm)	Mag. range
DFPL0.5×	0.050	198	0.32 - 2.85
DFPL0.75×	0.075	113	0.47 - 4.3
DFPLA01×-2	0.100	87.5	0.63 - 5.7
SZX-ACH1×	0.100	90	0.63 - 5.7
DFPL1.5×	0.150	53	0.95 - 8.55
DFPL2×-3	0.200	34	1.26 - 11.4

Eyepieces



Eyepiece	F.N.	Reticle diam.
WHS10×-H	22	Ø24
CROSS-WHS10×	22	-
WHS15×-H	16	Ø24
CROSS-WHS15×-H	16	-
WHS20×-H	12.5	Ø24
CROSS-WHS20×-H	12.5	-
WHS30×-H	7	Ø24

Advanced Ergonomic Design For Greater Working Comfort.

With elbows on the desk, operation is even easier.

When focus adjustment and zoom change can both be performed without removing the user's elbows from the desk, fatigue is significantly reduced. The transmitted light illuminator has been designed to be longer in front and back, but with reduced width and depth: this provides a generous 200 X 200mm area on which to place a large specimen, improves access to the specimen and also saves desk space.

Slimline body and 30-degree head tilt angle for better observation.

Because the zoom body is slim, there is better recognition of the specimen beneath the objective. The binocular head SZKBI30 and the trinocular head SZX-TR30 are both designed to be tilted at an angle of 30 degrees: the consequent improvement in working comfort, and the provision of adequate distance between the eyes and the specimen, combine to make observation easier.


Tilting binocular head and eyepoint adjuster for comfortable observation.

To make observation as comfortable as possible, Olympus has designed the tilting binocular head SZX-TBI to tilt from 5 to 45 degrees, and an eyepoint adjuster SZX-EPA. These allow different eyepoints. Users no longer have to bend their necks awkwardly to look into the microscope.



Built-in magnification display with integral click-stop.

The zoom handle on the right side of the zoom body is marked with a clearly visible magnification display (for use with the combination of 1 X objective and 10 X eyepiece). A click-stop mechanism (which can be overridden) easily gives the magnification sought. The exclusive magnification display rings allow easy reading of the magnification when objectives other than 1 X are being used.



ideal For Training And Group Discussions: Just Follow The Arrow.



Side by side discussion tube and arrow pointers are essential tools in classroom and training situations where aspects of operation and inspection must be communicated in detail. The arrow pointer can be displayed on images or monitors to highlight particular areas or features under discussion.

Side by side discussion tube / SZX-SD0

A generous 600mm of space separates primary and secondary observers, so there is no problem of crowding to compromise operational comfort. The integrated pointer can be displayed in different colors (yellow and red switchable) depending on the specimen.


SZX-SD0



SZX APT

Arrow pointer / SZX-APT

The arrow pointer can be superimposed on a picture or any observation image — a convenient, effective and economical system for classrooms and when making pictorial documents for quality assurance or quality control.



Triangular cantilever: red arrow (left), yellow arrow (right)

Additional Equipment For More Flexible, Efficient Image Recording.

When a trinocular head and a beam-splitter are combined a total of 3 photomicrography units or video cameras can be attached simultaneously. The SZX Series also features a range of mount adapters to accommodate to various additional image devices. Digital cameras are particularly useful, allowing easy image filing slide making and image transmission via a personal computer.

Beam splitter / SZX-BS

Combining this adapter and SZX-PHA or SZX-SLR allow a 35mm SLR camera, a video camera and another photomicrography unit to be attached on both sides of the body. The optical path is switchable in three stages: 100% observation, 100% photo/video, or 50% observation and 50% photo/video. The 100% optical path is best for observing and taking a dark specimen.



Photo adapter / SZX-PHA

Combining a wide range of camera mounts allows SZX microscopes to be linked to various video cameras.



SLR adapter / SZX-SLR

This unit allows a 35mm SLR camera to be attached. A separate OM mount adapter (PM-CAMS2) is necessary to attach the camera itself.



SZX-SLR+SZX-BS+SZX-PHA



SLX BS+SZX-PHA

A range of video camera adapters

A range of adapters is available to accommodate different video cameras to different intermediate magnifications and mounts. The 0.25 X C-mount video port U-TV0.25 X C, with rotating mount and

parfocal adjusting, combined with a 1/3" video camera, gives an image inscribed within the observation field. These adapters can be attached to the trinocular head SZX-TR30 and the photo adapter SZX-PHA directly.



Automatic Photomicrography system / PM30

With the PM30, shooting is simple — requiring nothing more than a touch on the shutter release — because the system's integrated super FL auto mode allows fully automatic exposure of fluorescent photos for consistently excellent results. Normal auto mode should be selected for automatic shooting of non-fluorescent specimen images. Among its wealth of features, the PM30 offers auto bracketing, 0.1% to 1% spot measurements, 3 types of photometry areas with a 30% average measurement, and data backup by memory card.



Automatic Photomicrography system / PM20

The PM20 system integrates FL auto mode for automatic exposure of fluorescent photomicrography with normal auto mode. More high-performance features are included, yet the design is actually reduced in size. Switching is available between 1% spot measurement and 30% average measurement, and up to 4 predetermined settings can be stored in the internal memory.



Automatic photography system / PM-10AK3

This allows the attachment of large-format cameras as well as 35mm cameras. The system is compact, yet packed with features that are easy to select and operate — including automatic exposure and auto/manual film advance.



PM30 specifications

Automatic exposure control unit	Measuring system	Two-dimensional split measuring
	Photometric modes	Super FL Auto mode, FL Auto mode, Auto mode, Manual mode, Time mode (AE lock, Multi Exp. also available)
	Auto exposure adjustment range	1/125sec. - 68min. (Auto mode, ISO100)
	Measuring area	0.1%, 1% spot measurements 30% average measurement
Acceptable cameras		35mm camera, 4" × 5" intermediate adapter, 3 1/4" × 4 1/4" Polaroid camera
Control unit	Data input via sheet switches and jog dial. Automatic ISO setting via DX code (with PM-C35DX mounted), manual setting also possible. Reciprocity failure adjustment function. Data backup by internal memory. Auto bracketing (3/5/7 frames).	
	Large LCD (backlit) screen (320 × 240 dot matrix) RS-232C interface IC memory card/printer interface (optional)	
Power supply		100 - 120V, 220 - 240V, 50/60Hz, 150VA

PM20 specifications

Automatic exposure control unit	Measuring system	Real-time measuring (TTL)
	Photometric modes	FL Auto mode, Auto mode, Manual mode, Time mode (AE lock, Multi Exp. also possible)
	Auto exposure adjustment range	1/125sec. - 68min. (Auto mode, ISO100)
	Measuring area	0.1% spot measurement, 30% average measurement
Acceptable cameras		35mm camera, 4" × 5" intermediate adapter, 3 1/4" × 4 1/4" Polaroid camera
Control unit	Data input via sheet switches and jog dial. Automatic ISO setting via DX code (with PM-C35DX mounted), manual setting also possible. Reciprocity failure adjustment function. Data backup by internal memory Auto bracketing (3/5/7 frames)	
	LCD (backlit) screen	
Power supply		100 - 120V, 220 - 240V, 50/60Hz, 50VA

PM-10AK3 specifications

Automatic exposure control unit	Measuring system	Real-time (TTL)
	Exposure modes	AUTO and TIME modes
	Auto exposure adjustment range	1/125 sec - 4 minutes (Auto mode, ISO100)
	Measuring area	30% average measurement
Acceptable cameras		35mm camera, 3 1/4" × 4 1/4" Polaroid camera, 4" × 5" sheet film (Recordata back can be mounted on 35mm camera back)
Control unit		Data input via dial, manual ISO setting, estimated exposure time is indicated in 4 step LED

Seen In The Right Light: Illumination That Suits The Specimen.

■ **Reflected Light Illuminators** (Different types of light sources and light guides may be offered in each area)

Light guide power supply /

LG-PS2

This long-life 12V 1 OOW halogen lamp assures bright, reliable illumination for long hours, and there's a mechanical stop which allows light intensity adjustments without changing the color temperature. The unit is small and compact, easy to carry, and conforms to relevant safety standards, making it the ideal light source for the SZX series.

Coaxial illuminator /

szx-ILLC

Used with the dual flexible light guide LG-DF, this illuminator provides bright, steady illumination without tiresome centering adjustments. It's particularly effective for observing structure, such as imperfections on metal surfaces, patterns on IC or LCD, which don't show up under conventional oblique illumination.

Dual combination light guide /

LG-DFI

The standard SZX light guide can be mounted directly on SZX-FOF/SZX-FO and keeps the observation position properly illuminated even if focusing shifts when one specimen is changed for another. Mounting is designed to avoid hindering observation.

Dual inter-lock light guide /

LG-DI

Employs random positioning to give bright, steady illumination — especially effective when high contrast images are required. The condenser lens HLL301 can be mounted on the end of the light guides.



LG-DI + LG-PS2

Ring light guide /

LG-R66

With its 066 mounting diameter, this ring light guide has been specially developed for stereo microscope compatibility. When mounted with ring light adapter SZX-LGR66, it provides bright, clear images by illuminating from several angles, thus eliminating obstructive specimen shadows.

Reflected illumination adapter /

LG-ILLR (available from March 1998)

This adapter is attached to ring light guide LG-R66. This ensures soft, effective illumination even when conducting observation of highly reflective specimens.



SZX12 + SLX-ILLC + LG DF + LG PS2



SZX12 + LG-DFI + LG-PS2



SZX12 + SLX-LGR-66 + LGR66 + LG-PS2

■ Transmitted Light Illuminators

Transmitted light illuminator / SZX-ILLK

This cost-effective transmitted light illuminator with built-in 6V 30W halogen lamp is particularly suitable for use with 1X objectives. It uses oblique incident light illumination to provide contrasting images of a transparent specimen.



High-level transmitted light illuminator / SZX-ILLB

With its variable magnification range (from 1X to high magnification), this unit provides clear, effective illumination. Light volume and color temperature are adjusted by means of three built-in filters. Used with a high magnification objective such as the SZX-AL20X, it permits observation of extremely detailed, high contrast structures under unique, slit stop oblique illumination — an advantage which cannot be matched by conventional stereo microscopes.



Brightfield/darkfield transmitted illuminator / SZX-ILLD

This unit features newly developed optical characteristics for bright illumination distributed evenly over a wide field diameter. The unit can switch between brightfield illumination optics which illuminate every area in the lowest magnification field diameter (Ø63mm) in an DFPLFL0.5 X PF objective, to darkfield illumination optics for a 5 X total magnification field diameter (Ø44mm). Wide range of general applications.




Transmitted light guide adapter / SZX-TLGAD Single flexible light guide/LG-SF*

This combination allows the light source to be separated from the transmitted illuminator, so that there is no temperature increase on the stand at all. Indispensable in cases when the specimen temperature must be strictly managed.




* Different types of light guides may be offered in each area.



Pattern on the Charge Modulation Device (CMD)



Surface of the rack after resistance test



Diode on the circuit board

A Wide Range Of Accessories For Other Observation Applications.



Large stand / SZX-STL

The stand offers excellent stability, making it ideal for photomicrography and video documentation.




Universal stand / SZ-STU2

This stand is designed for observation and photomicrography of specimens which are too large for a stand. A structure based on dual horizontal poles and linear ball bearings ensures smooth horizontal movement and rotation. The mechanism can be tilted forward, backward, right and left, allowing quick and precise approach to the chosen observation area.




BX stage adapter type 1 / SZX-STAD1

Allows use of the BX rotating stage (U-SRG) to various SZX stands and transmitted light illuminators. This is especially valuable in polarized observations and photomicrographic framing.



BX stage adapter type 2 / SZX-STAD2

Allows use of a BX mechanical stage by combining with various SZX stands and transmitted light illuminators. Particularly suitable for accurate X-Y movement of specimens.



BH stage adapter type 1 / SZH-STAD1

Allows use of a BH2 mechanical stage (BH2-SH) by combining with various SZX stands and transmitted light illuminators. Particularly suitable for accurate X-Y movement of specimens.




Fluorescence units



Drawing attachment / SZX-DA

This attachment enables the user to make an accurate drawing of the specimen — a valuable alternative to photomicrography, since only the required areas are drawn. Can be mounted on either the right or the left side.



Debris on the wafer

Fluorescence units

(available from March 1998)

The fluorescence units make it possible to inspect subsurface imperfections or to observe photoresist residues under fluorescence microscopy.



Simple polarizer/SZX-PO

The polarizer is mounted on the stage plate mount of the transmitted light illuminator, and used to observe complex images such as thin mineral substances. During observation, an optional analyzer (see below) should be mounted on the objectives.



Large stage plate/SZX-CL

(available from March 1998)

The large plate can be attached with 200 X 200mm glass plate SZX-CL to the transmitted light unit. The plate covers almost the whole surface of the stand and makes it easy to clean or sterilize the stage. It's also possible to insert a plane illuminator, or sheets of colored paper, to obtain a variety of lighting effects.



Focusing unit/SZ-FO

This unit is movable, allowing simultaneous focus on itself and on the specimen. As a result, the user can observe or photograph specimens of different thickness without changing the eyepoint.



Rotatable analyzer / SZX-AN


This analyzer is mounted on the top of the objectives and used in conjunction with the simple polarizer SZX-PO. Used for both polarized transmitted light and polarized reflected light observation, the analyzer is also effective in avoiding direct reflected light.

A Full Range Of Effective Combinations To Suit Every Need

Examples of SZX9 combination

SZX9-1213

This economical combination features 30 degree binocular head SZX-BI30, focusing unit SZX-FO and apochromat objective SZX-ACH1 \times , and is suitable for a wide variety of basic observations and work applications. The SZX-ACH1 \times objective has a 90mm working distance.




Lead frame



SZX9-2212


Combination with tilting binocular head SZX-TBI, focusing unit SZX-FO and apochromat objective DFPLAPO1 \times -2. The tilting tube, which enables to change the tilt angle of eyepieces from 5 to 45 degrees, allows extended observations without fatigue, while the apochromat objective ensures natural color reproduction.



Wire bonding section on the AF sensor

SZX9-3112

This combination consists of trinocular tube SZX-TR30, coarse and fine focusing unit SZX-FOF, apochromat objective DFPLAPO1 \times -2 and arrow pointer SZX-APT. The SZX-FOF unit enhances operational efficiency by making delicate focusing easy, while the arrow pointer (compatible with photo, monitor and video printer), helps to ensure efficient image recording.



Expanded surface of the rack



SZX-1 2 1 3

■ Observation head, intermediate tube

- ① SZX-BI30 ② SZX-TBI
- ③ SZX-TR30 ④ SZX-BI30+SZX-BS
- ⑤ SZX-TBI+SZX-BS
- ⑥ SZX-BI30+SZX-DO

■ Focusing unit

- ① SZX-FOF
- ② SZX-FO

■ Objective

- ① DFPLAPO1×PF ② DFPLAPO1×2 ③ SZX-ACH1×
- ④ DFPLAPO1×PF+DFPLFLO0.5×PF
- ⑤ DFPLAPO1×PF+DFPLFLO1.6×PF
- ⑥ DFPLAPO1×PF+DFPLAPO1.6×PF+SZX-AL20
- ⑦ DFPLAPO1.2×PF ⑧ DFPLAPO1×PF+SZX-AL20×


■ Stand

- ① SZX-ST ② SZX-ILLK ③ SZX-ILLB ④ SZX-ILLD
- ⑤ SZX-ILLC ⑥ SZX-STL2(SDO) ⑦ SZX-STL2(SDO)

Examples of SZX12 combination

SZX12-2111

A package for standard applications: tilting binocular head SZX-TBI, coarse and fine focusing unit SZX-FOF and apochromat parfocal objective DFPLAPO1×PF. This combination designed for fatigue-free observation via the tilting binocular head. Together with an image recording intermediate tube such as the SZX-BS, this set can be connected to photomicrography and video units without using a trinocular head.




Pattern on the Charge Modulation Device (CMD)



SZX12-3111

This set consists of 30 degree trinocular head SZX-TR30, coarse and fine focusing unit SZX-FOF and apochromat parfocal objective DFPLAPO1×PF. The trinocular head lets you use up to three photomicrography and video units when combined with an image recording intermediate tube like the SZX-BS.




Pattern on the Charge Modulation Device (CMD)



SZX12-3115

This versatile combination features 30 degree trinocular head SZX-TR30, coarse and fine focusing unit SZX-FOF, apochromat parfocal objective DFPLAPO1×PF, objective DFPLAPO1.6×, revolving nosepiece SZX-2RE and arrow pointer SZX-APT. The revolving nosepiece and two objectives (other combinations are also possible) allow a wide range of magnification changes. The trinocular head and arrow pointer make a highly functional combination, especially in terms of image recording; in addition, the arrow pointer can be imprinted on to film for more effective presentations and training sessions.



Wire bonding section on the AF sensor



Specifications

Item	Specifications					
	SZX-ZB12			SZX-ZB9		
Zoom microscope bodies	Zoom variable magnification system with parallel optical axis, Zoom drive system: Horizontal handle Click-stop for various zoom positions incorporated					
Total magnification (Objective 1×, eyepiece 10×)	Zoom ratio: 12.86(0.7×-9×) Total magnification indication: 7/10/12.5/16/20/25/32/40/50/63/90			Zoom ratio: 9.0 (0.63×-5.7×) Total magnification indication*: 6.3/8/10/12.5/16/20/25/32/40/50/57		
	Objective mounting: Screw mount					
	Built-in aperture iris diaphragm			Aperture diaphragm unit (SZX-AS) is available		
Focusing assembly	SZX-FOF: Fine focusing unit Focus: Rack and pinion with roller guide(with torque adjustment ring for coarse focusing), built-in counter balance, coarse and fine coaxial handle, coarse handle stroke: 80mm (stroke per rotation: 36.8mm), fine handle stroke: 80mm (stroke per rotation: 1.5mm)					
	SZX-FO: Focusing unit Focus: Rack and pinion with roller guide(with torque adjustment ring for coarse focusing), counter balance optional, coarse handle, coarse handle stroke: 80mm, coarse handle stroke per rotation: 21.2mm					
Observation heads	SZX-TR30: 30 degree trinocular head Tilting angle: 30°, Light path selection: 2 steps (100% binocular tube, 20% binocular tube and 80% photo), Interpupillary distance adjustment: 50 - 76mm, with eyepiece fixing knob, eyepiece: WHS series					
	SZX-BI30: 30 degree binocular head Tilting angle: 30°, Interpupillary distance adjustment: 50 - 76mm, with eyepiece fixing knob, eyepiece: WHS series					
	SZX-TBI: Tilting binocular head Tilting angle: 5° - 45°, Interpupillary distance adjustment: 50 - 76mm, with eyepiece fixing knob, eyepiece: WHS series					
Stands	SZX-ST: Stand Pillar height: 270mm, base dimension: 300(W)×260(D)×30(H)mm, stage clips are mountable, with stage adapter fixing screw holes					
	SZX-STL: Large stand Pillar height: 400mm, base dimension: 400(W)×350(D)×28(H)mm, stage clips are mountable, drop prevention collar (SZX-R) is necessary					
Objectives W.D.: Working Distance	For SZX-ZB12			For SZX-ZB9		
	Objectives	N.A.	W.D.	Objectives	N.A.	W.D.
PF: Parfocal objectives	DFPLFL0.3×	0.036	130mm*	DFPL0.5×	0.05	198mm*
Auxiliary long pillar and drop prevention collar (SZX-R) are necessary	DFPLFL0.45×	0.05	198mm	DFPL0.75×	0.075	113mm
	DFPLFL0.5×PF	0.055	70mm	DFPLAPO1×-2	0.10	87.5mm
	DFPLAPO1×PF	0.110	74mm	SZX-ACH1×	0.10	90mm
	DFPLAPO1.2×PF	0.13	60mm	DFPL1.5×	0.15	53mm
	DFPLFL1.6×PF	0.176	34mm	DFPL2×-3	0.2	34mm
	SZX-AL20×	0.275	107mm			
Eyeiece *Possible to insert micrometer (ø24mm, 1.5mm thick)	WHS10×-H * F.N. 22 CROSS-WHS10× F.N. 22, with cross lines WHS15×-H * F.N. 16 CROSS-WHS15× F.N. 16, with cross lines WHS20×-H * F.N. 12.5 CROSS-WHS20× F.N. 12.5, with cross lines WHS30×-H * F.N. 7					

Reflected light illuminators specifications (Different types of light sources and light guide may offered in each area)

	SZX-ILLC	LG-DF/DI	LG-R66
Specifications	Magnification factor: 1.5× Light source: LG-PS2 Light guide: LG-DF Fiber diameter: ø6 Length: 1,000mm Simple polarized light observation is available	Dual inter-lock light guide Light source: LG-PS2 Fiber diameter: ø8 LG-DF: Flexible part 900mm, Inter-lock part 500mm LG-DI: Inter-lock part 500mm	Ring light guide Light source: LG-PS2 Fiber diameter: ø8 Flexible part 1,000mm Fixing diameter 66mm, Minimum W.D. 30mm Ring light guide adapter for SZX: SZX-LGR66
Option	—	HILL301: spot lens LG-FAD: ø25 filter adapter (worldwide releases March 1998)	LG-R66PL: Polarizer/analyzer set for LG-R66 LG-ILLR: Reflected illumination adapter (worldwide releases March 1998)
Features	Brighter and evener illumination without adjustment	LG-DF: Fixed to the focusing unit and illumination is synchronizing with the movement of the microscope.	Bright, even, and no shading illumination is achieved in case of any objectives.
Light source specifications	Lamp: 12V 100W long-life (1,000 hrs.) halogen lamp with reflected mirror Rated voltage: 100-120V/220V-240V 50/60Hz Power consumption: 350VA Protection functions: Overheat detection by temperature sensor, over current detection of input/output current Other functions: Light intensity control and lamp ON/OFF control by external signal, detection of burned out lamp Dimensions: 120(H) × 120 (W) × 235(D)mm Weight: approx. 1.5kg Option: Color filter		

Transmitted illuminators specifications

Item	Specifications		
	SZX-ILLK	SZX-ILLB	SZX-ILLD
Light source	PHILLIPS 5761 6V30W/40V halogen lamp (average lamp service life: approx. 100 hours under rated usage conditions)		
Light intensity adjustment	Continuously variable system (with built-in transformer)		
Effective illuminated area	ø40mm		Brightfield: ø63mm, Darkfield: ø45mm
Built-in filter	—	LBD, ND6, ND25 one for each	FR, LBD, ND25 one for each
Add-on filter	ø45mm frosted filter (45OPL), provided		
Illumination mode	Transmitted brightfield illumination Oblique illumination	Transmitted brightfield illumination Oblique illumination, Slit aperture illumination	Transmitted brightfield illumination Transmitted darkfield illumination
Magnification selection	—	2-step selection between 1× - 1.2× and 1.5× or higher objectives	—
Pillar height	270mm		
Weight	Approx. 6kg (13.2lb)	Approx. 6.2kg (13.6lb)	
Rated voltage	100V area: 100/110-120V 0.6A 50/60Hz 200V area: 220/230-240V 0.3A 50/60Hz		

Total magnifications and actual field diameters of SZX-ZB12


Objective	Eyepiece							
	WHS10×-H		WHS15×-H		WHS20×-H		WHS30×-H	
	total mag.	field diameter	total mag.	field diameter	total mag.	field diameter	total mag.	field diameter
DFPLFL0.3×	2.1× - 27×	104.8 - 8.1mm	3.15× - 40.5×	76.2 - 5.9mm	4.2× - 54×	59.5 - 4.6mm	6.3× - 81×	33.3 - 2.6mm
DFPLFL0.45×	3.15× - 40.5×	69.8 - 5.5mm	4.73× - 60.7×	50.8 - 4.0mm	6.3× - 81×	39.7 - 3.1mm	9.6× - 121.5×	22.2 - 1.7mm
DFPLFL0.5×PF	3.5× - 45×	62.9 - 4.9mm	5.25× - 67.5×	45.7 - 3.6mm	7× - 90×	35.7 - 2.8mm	10.5× - 135×	20.0 - 1.6mm
DFPLAPO1×PF	7× - 90×	31.4 - 2.4mm	10.5× - 135×	22.9 - 1.8mm	14× - 180×	17.9 - 1.4mm	21× - 270×	10 - 0.78mm
DFPLAPO1.2×PF	8.4× - 108×	26.2 - 2.0mm	12.6× - 162×	19.0 - 1.5mm	16.8× - 216×	14.9 - 1.2mm	25.2× - 324×	8.3 - 0.65mm
DFPLFL1.6×PF	11× - 144×	20 - 1.5mm	16.8× - 216×	14.3 - 1.1mm	22.4× - 288×	11.2 - 0.87mm	33.6× - 432×	6.25 - 0.49mm
SZX-AL20×*	100× - 225×	22 - 0.98mm	150× - 337.5×	1.6 - 0.71mm	200× - 450×	1.25 - 0.56mm	300× - 675×	0.44 - 0.31mm

* This is the auxiliary objective mounted on DFPLAPO1×PF. Figures indicate recommended magnifications.


Total magnifications and actual field diameters of SZX-ZB9

Objective	Eyepiece							
	WHS10×-H		WHS15×-H		WHS20×-H		WHS30×-H	
	total mag.	field diameter	total mag.	field diameter	total mag.	field diameter	total mag.	field diameter
DFPL0.5×	3.15× - 28.5×	69.8 - 7.7mm	4.7× - 42.8×	50.8 - 5.6mm	6.3× - 57×	39.7 - 4.4mm	9.5× - 85.5×	22.2 - 2.5mm
DFPL0.75×	4.7× - 43×	46.6 - 5.1mm	7.1× - 64.1×	33.9 - 3.7mm	9.5× - 85.5×	26.5 - 2.9mm	14.2× - 128×	14.8 - 1.6mm
DFPLAPO1×-2	6.3× - 57×	34.9 - 3.9mm	9.5× - 85.5×	25.4 - 2.8mm	12.6× - 114×	19.8 - 2.2mm	18.9× - 171×	11.1 - 1.2mm
SZX-ACH1×	6.3× - 57×	34.9 - 3.9mm	9.5× - 85.5×	25.4 - 2.8mm	12.6× - 114×	19.8 - 2.2mm	18.9× - 171×	11.1 - 1.2mm
DFPL1.5×	9.5× - 85.5×	23.3 - 2.6mm	14.2× - 128.3×	16.9 - 1.9mm	19× - 171×	13.2 - 1.5mm	28.4× - 256.5×	7.4 - 0.82mm
DFPL2×-3	12.6× - 114×	17.5 - 1.9mm	18.9× - 171×	12.7 - 1.4mm	25.2× - 228×	9.9 - 1.1mm	37.8× - 342×	5.6 - 0.61mm


SZX12 + SZX-ST dimensions




SZX9 + SZX-ST dimensions




SZX12 + SZX-ILLC dimensions




SZX9 + SZX-ILLC dimensions




SZX12 + SZX-ILLD dimensions




SZX9 + SZX-ILLD dimensions




Photomicrography combination dimensions (with SZX-BS)




Photomicrography combination dimensions (SZX-TR30 + 35mm camera back)




SZX12 + SZX-STL dimensions




SZ-STU2 dimensions



SZX12 + SZX-SDO dimensions




SZX12 + SZX-DA dimensions




Unit: mm

SZX SYSTEM DIAGRAM



*1 To be released in March, 1990 *2 Exclusive for the SZX-ZB9 *3 Incorporated in SZX-ILLC *4 Cannot be attached with SZX-TR30 *5 Cannot be used with SZX-SDO *6 Different types of light



Accessories and light guides may be offered in each area. *7 SZH-P400 is the standard feature of SZX-STL and SZX-STL2.

Web site addresses: <http://www.olympus.co.jp>
<http://www.olympus.com>
<http://www.olympus-europa.com>



Specifications are subject to change without any obligation on the part of the manufacturer

Olympus business areas

	Medical and health-care area
	Imaging and information area
	Industrial applications area



OLYMPUS OPTICAL CO., LTD.
 2-43-2 Hongo-ya Bunkyo-ku Tokyo Japan
 OLYMPUS OPTICAL CO. (EUROPA) GMBH.
 Postfach 10 49 08 20024 Hamburg Germany
 OLYMPUS AMERICA INC.
 3 Corporate Center Drive Melville NY 11747-3157 U.S.A
 OLYMPUS SINGAPORE PTE. LTD.
 111 Henderson Road #13 03 Henderson Industrial Park Singapore 159552 Singapore
 OLYMPUS OPTICAL CO. (UK) LTD.
 28 Harsnett Street London EC1A 3DT United Kingdom
 OLYMPUS AUSTRALIA PTY. LTD.
 104 Fernside Gully Road Oakleigh VIC 3166 Australia